

superlite



FEATURES

- Hardwired to high resolution L7HD probe
- DuraScan® case: water, dust and impact resistant
- 3.3 lb (1.5 kg)
- B, B/M modes
- Auto enhanced field of view (auto adjusts focus, frequency, gain and other parameters)
- Digital beam-former
- Gestational tables
- Stores still frames and video loops
- 4GB storage for over 10,000+ still frames
- Voice label your stored images
- Li-Ion smart batteries with 6+ hours run time
- Multiple headset options
- 3 year warranty (warranty extension available)

DESCRIPTION

With 90% of the features and the same image quality of the Ibex® Pro and Pro/r, the Ibex® SuperLite is the ideal choice for bovine ultrasound. Hardwired to the L7HD, our latest high resolution transducer, the SuperLite is designed specifically for bovine reproduction, both dairy and beef.

The Ibex® SuperLite has all the features you need in an ultra light, ruggedized package. The new L7HD probe provides high definition images usually found in higher end ultrasound units.

With most of the features of the Ibex® Pro and Pro/r but without the added weight and cost of the screen, the Ibex® SuperLite is ideal for bovine practice and outdoor field work. Video headsets are great for working in both bright sunlight and low light conditions—worn around the waist, the unit is lightweight and makes for streamlined workflow. Like the other products in the Ibex® line, the SuperLite is rugged and water-resistant, with a long battery life.

