

superlite/c



FEATURES

- Hardwired to high resolution CLi4HD probe
- DuraScan® case: water, dust and impact resistant
- 3.3 lb (1.5 kg)
- B, B/M modes
- Auto enhanced field of view (auto adjusts focus, frequency, gain and other parameters)
- Digital beam-former
- Gestational tables
- Stores still frames and video loops
- 4GB storage for over 10,000+ still frames
- Li-Ion smart batteries with 6+ hours run time
- Multiple headset options
- 3-year warranty (warranty extension available)

DESCRIPTION

The IBEX® SuperLite/c is the ideal choice for large animal ultrasound. Hardwired to the CLi4HD, our latest high resolution transducer, the SuperLite/c is designed specifically for bovine reproduction, both dairy and beef.

The IBEX® SuperLite/c has all the features you need in an ultra light, ruggedized package. The new CLi4HD probe provides high definition images usually found in higher end ultrasound units.

With most of the features of the IBEX® PRO series but without the added weight and cost of the screen, the IBEX® SuperLite/c is ideal for bovine practice and outdoor field work. Video headsets are great for working in both bright sunlight and low light conditions—worn around the waist, the unit is lightweight and makes for streamlined workflow. Like the other products in the IBEX® line, the SuperLite/c is rugged and water-resistant, with a long battery life.

