

IBEX® PRO/c

Designed specifically for bovine and equine reproduction. New technology provides clear, high definition scans seen only on higher end ultrasound units.



Features

- Hardwired to high definition **CLi4HD** probe for improved resolution and ergonomics
- DuraScan® case: water, dust and impact resistant
- 5.5 lb (2.5 kg)
- B, B/B, B/M and M scan modes
- Auto enhanced field of view (auto adjusts focus, frequency, gain and other parameters)
- Digital beam-former
- Gestational table
- Stores still frames and video loops
- 4GB storage for over 10,000+ still frames
- Li-Ion smart batteries with 4+ hours of runtime
- Multiple video headset options available
- 3-year warranty (warranty extension available)

Features

- Hardwired to high resolution **CLi4HD** probe for improved resolution and ergonomics
- DuraScan® case: water, dust and impact resistant
- 3.3 lb (1.5 kg)
- B, B/M modes
- Auto enhanced field of view (auto adjusts focus, frequency, gain and other parameters)
- Digital beam-former
- Gestational tables
- Stores still frames and video loops
- 4GB storage for over 10,000+ still frames
- Li-Ion smart batteries with 6+ hours of run time
- Multiple video headset options
- 3-year warranty (warranty extension available)

IBEX® SuperLite/c

The LITE/c has all the features you need in an ultra light, ruggedized package. The new hardwired **CLi4HD** probe provides high definition images usually found in higher end ultrasound units. Enhanced battery performance providing 6+ hours of continuous runtime.

