



E.I. Medical Imaging®

New products, continuous innovation

IBEX® PRO/r

Designed specifically for bovine and equine reproduction. New technology provides clear, high definition scans seen only on higher end ultrasound units.



Features

- Hardwired to high definition L7HD probe—improved resolution and ergonomics
- DuraScan® case: water, dust and impact resistant
- 5.5 lb (2.5 kg)
- B, B/B, B/M and M scan modes
- Auto enhanced field of view (auto adjusts focus, frequency, gain and other parameters)
- Digital beam-former
- Gestational table
- Stores still frames and video loops
- 4GB storage for over 10,000+ still frames
- File management system
- Li-Ion smart batteries with 4+ hours of runtime
- Multiple headset options available
- 3 year warranty (warranty extension available)

Features

- Hardwired to high resolution L7HD probe—improved resolution and ergonomics
- DuraScan® case: water, dust and impact resistant
- 3.3 lb (1.5 kg)
- B, B/M modes
- Auto enhanced field of view (auto adjusts focus, frequency, gain and other parameters)
- Digital beam-former
- Gestational tables
- Stores still frames and video loops
- 4GB storage for over 10,000+ still frames
- Voice label your stored images
- Li-Ion smart batteries with 6+ hours of run time
- Multiple headset options
- 3 year warranty (warranty extension available)

IBEX® SuperLite

The Ibex® SuperLite has all the features you need in an ultra light, ruggedized package. The new hardwired L7HD probe provides high definition images usually found in higher end ultrasound units. Enhanced battery performance providing 6+ hours of continuous runtime.



E.I. Medical Imaging • 866.365.6596 • eimedical.com

