

Next-Generation Probes – Wired or Wireless Freedom

Designed for Dog Breeders

See Pregnancies Clearly. Scan with Confidence.

The **ScanX® Portable Veterinary Ultrasound System** is a lightweight, easy-to-use ultrasound designed specifically for **small animal pregnancy detection**, with dog breeders in mind. Now upgraded for 2026, ScanX introduces **next-generation probes** that allow you to scan **wired or wirelessly**, giving you total flexibility in the kennel, whelping room, or on the go.

Safe, reliable, and beginner-friendly, ScanX delivers **clear, high-resolution images** so you can confidently confirm pregnancies without stress.

What's New for 2026

Detachable Wired + Wireless Probe Technology

The new **ScanX® next-generation probe** gives you the best of both worlds:

- **Scan wirelessly** for complete freedom of movement
- **Attach the cable** when you prefer a traditional, secure connection
- **Same image quality** in both wired and wireless modes
- Lightweight, balanced, and comfortable to hold

No compromises. Just scan the way that works best for you.

Perfect for Dog Breeders

- Designed for **canine pregnancy detection**
 - Simple touchscreen operation
 - Easy to learn, even for first-time ultrasound users
 - Quiet, portable, and kennel-friendly
 - Used by breeders, rescue centres, and professionals worldwide
-

What's Included

- **Microsoft Surface Go 2 Tablet (10.5" touchscreen)**
 - **ScanX Ultrasound Software** with obstetrics presets
 - **One Ultrasound Probe** (Micro-Convex or Convex)
 - **Automatic cloud backup** of images and videos
 - North American service, support, and remote assistance
-

Choose the Right Probe for Your Dogs...

Micro-Convex Probe (5–8.5 MHz)

Next-Generation Detachable + Wireless (Standard)

Best for **small to medium dog breeds** and leaner dogs.

Ideal for:

Chihuahua, Toy & Standard Poodle, Corgi, Cocker Spaniel, Labrador, German Shepherd, Border Collie, Doberman, Dachshund, Yorkshire Terrier, Whippet, Shih Tzu, and similar breeds.

- Small footprint for easy positioning
 - Higher frequency for excellent detail
 - Scanning depth up to 15 cm
-

Convex Probe (2.5–5 MHz)

Traditional Wired Probe

Best for **large and giant dog breeds** with more body mass.

Ideal for:

Saint Bernard, Newfoundland, Great Dane, Bernese Mountain Dog, Old English Sheepdog, and similar large breeds.

- Larger footprint
 - Deeper penetration
 - Scanning depth up to 30 cm
-

Key Features at a Glance

- 10.5” high-resolution touchscreen
 - Lightweight system (approx. 1.2 lbs)
 - Battery life up to **10 hours**
 - B, M, 2B, B/M, and 4B scan modes
 - Built-in Wi-Fi
 - USB ports for exporting images and videos
 - Compatible with Microsoft Surface accessories
 - Automatic cloud image & video backup
-

Reliable Support You Can Trust

If you ever need help, you’re never on your own.

ScanX can be **serviced remotely**, and future software updates are completed **without sending your unit away**

Scan Smarter. Scan Easier. Scan Your Way. - ScanX® Portable Ultrasound – Built for Breeders