ViiSona
ultrasound & education

PORTABLE VETERINARY ULTRASOUND

IBEX® PRONXT BY E.I. MEDICAL IMAGING

ProNXT

High-definition bovine and equine repro ultrasound with HerdIQ™ AI

Maximum image detail for bovine and equine reproduction: a 25% larger sunlight-readable display, up to 60 fps, full measurement menus and HerdIQ™ AI on board.

- **25% larger sunlight-readable display**
- **Up to 60 fps live imaging**
- **HerdIQ™ AI on board**



STARTING AT

ProNXT with eCLi6 probe — **$20,900**
Display and streaming built in

InSite NXT OLED headset — **$3,040**
Optional

60 fps
Up to, for a smooth live image

25%
Larger high-resolution display

4 h
Continuous scanning per battery

eCLi6
Hybrid probe, 3–9.5 MHz

www.viisona.com · 1-800-672-9182

Ideal for

Cattle and horses.

ADDITIONAL APPLICATIONS

Cow/calf, dairy and equine practices that want maximum image detail.

01	Bovine and equine reproduction
02	Ovaries, uterus and early pregnancy
03	Pregnancy diagnosis with HerdIQ AI
04	Fetal aging with HerdIQ Aging AI
05	Lung exams
06	Detailed measurements

The next generation of large-animal repro

Bigger, brighter display

A 25% larger, high-resolution, sunlight-readable display for easy viewing at the chute or in the stall.

Faster live image

Up to 60 frames per second for a smooth, responsive picture as you scan.

eCLi6 hybrid probe

Linear and curved linear in one hardwired probe: the resolution of a flat linear with the depth and wider view of a convex, arms-in or arms-free.

Redesigned image software

Tissue harmonics, synthetic imaging and speckle reduction for a significant jump in image quality.

Simple, sealed controls

A simplified sealed keyboard and touchpad, easy to clean, with full measurement menus.

Stream and scan all day

Built-in IBEXStream to phones and tablets. 4 hours per battery, auto stand-by saves 20%, field-swappable with an external dock.

DIAGNOSTIC IMAGES



ProNXT with hardwired eCLi6 probe



Bovine pregnancy image



Simplified sealed keyboard and touchpad

HERDIQ™ AI · BUILT INTO YOUR NXT

A second set of eyes on every cow

Trained on tens of thousands of bovine ultrasound images, HerdIQ runs directly on the NXT and reads the image as you scan. No extra steps, no field connectivity, no change to your technique.

HerdIQ Pregnancy AI

- Shows "Pregnant" in real time as you scan
- 99% accurate from 28 days of gestation
- Lights up when it finds the uterus and ovaries
- Ideal for year-round dairy pregnancy checks

HerdIQ Aging AI · NEW

- Estimates gestational age live, with no measuring
- Typically within 3 days from 30 to 130 days; within 5 days to 150
- Green, yellow and red confidence lights
- Repro Workflow: log cows by visual or EID tag, export reports

HerdIQ supports your expertise; it doesn't replace it. A training partner for new vets and staff, with instant feedback on image quality.

HERDIQ PLANS
Beef: 3-month trial $870 or season package $1,570 · Dairy: $700/yr

Build your configuration

ITEM		STARTING PRICE
ProNXT with hardwired eCLi6 probe	Screen and IBEXStream included	**$20,900**
InSite NXT OLED headset (optional)	3× resolution, 70% larger image	**$3,040**
HerdIQ™ AI plans	Pregnancy and Aging AI	**From $700/yr**

Hardwired probe

eCLi6 hybrid curved linear (hardwired) — **3–9.5 MHz**

Resolution of a flat linear with the deeper reach (to 24 cm) and wider view of a convex. Bovine and equine reproduction and lung exams, arms-in or arms-free.

> *"Image quality is cutting edge, keyboard is very simple and intuitive, all at a price that works for cattle veterinarians."*
>
> **TK Livestock Services**

The set includes

- ProNXT with hardwired eCLi6 probe
- Field-swappable lithium-ion battery
- AC power supply
- Custom hard case

Optional accessories

- InSite NXT OLED headset
- Spare battery
- External charger dock (AC or car)
- I.C.E.® arms-free extension arm



Technical data

Display	25% larger, high-resolution, sunlight-readable
Frame rate	Up to 60 fps
Probe	Hardwired eCLi6 hybrid curved linear, 3–9.5 MHz, max depth 24 cm
Image processing	Tissue harmonics, synthetic imaging, speckle reduction
Controls	Simplified sealed keyboard and touchpad; full measurement menus
Streaming	Built-in Wi-Fi via IBEXStream, iOS and Android
Headset	Compatible with InSite NXT OLED headset (optional)
AI	HerdIQ™ Pregnancy and Aging AI (subscription), runs on the unit
Battery	Field-swappable, 4 h continuous; auto stand-by saves 20%; external charger dock, AC or car
Construction	IBEX rugged build, RoHS compliant
Cables	Kevlar-reinforced, tested to 1.2 million flex cycles
Software	Firmware updates for the life of the system
Made in	Designed, manufactured and supported in Colorado, USA

Technical service & support in Canada

Authorized service technicians

Our authorized service technicians maintain and repair your IBEX equipment here in Canada, with full technical service from our Montreal-area office.

Loaner systems

If your scanner needs service, we can lend you a replacement system so your practice keeps working without interruption.

New trade-in program

Trade in your current ultrasound system and trade up to a new IBEX system. Ask us for a trade-in evaluation.

Trade in. Trade up.

Ask us what your current system is worth toward an upgrade.

Flexible financing

Lease and financing options are available. Ask us for terms that fit your cash flow.

Manufacturer warranty

1-year manufacturer warranty · extended warranties available

Sample scans & info

Quotes, demos & service

ViiSona
ultrasound & education

SALES & TECHNICAL SERVICES
38 Rue Canvin
Kirkland, QC H9H 4S4
Canada

www.viisona.com · info@viisona.com · 1-800-672-9182

Starting prices in Canadian dollars (CAD), before applicable taxes and shipping. Final price depends on the probes and accessories selected. Specifications supplied by E.I. Medical Imaging and subject to change without notice. Pricing effective November 2026. Manufactured by E.I. Medical Imaging, Colorado, USA · Sold and serviced in Canada by ViiSona.